



Pregnancy test instructions

This is a diagram of your pregnancy test




- Remove all packaging
- Remove the cap from the test
- As you urinate, hold the absorbent tip of the test in the urine stream or collect some urine in a clean pot and dip the tip into the urine for 5 to 10 seconds



or



- Wait 5 to 10 minutes
- Read the pregnancy test



✓ This is **negative**.



✗ This is **positive**.
The abortion may not have worked.
Call BPAS immediately.



✗ This is **invalid**.
The abortion may not have worked.
Call BPAS immediately.

